RESTECH
As simple as breathing

RESMON PRO FULL SOFTWARE VERSION 23.1.0a update package for systems currently running version 23.0.0a.

UPDATE INSTRUCTIONS

Version 0 – 26th Aug 2026

IMPORTANT: This document details the steps to update a RESMON PRO FULL with software version 23.0.0a to version 23.1.0a. If the current installed software is older than 23.0.0a, these instructions are not applicable.

PREREQUISITES:

- USB Drive: You will need a USB stick with at least 200MB of free space.
- Format: Ensure the USB stick is formatted as FAT32.
- Power: Ensure the Resmon PRO FULL is connected to AC power throughout the entire process.

DOWNLOAD THE UPDATE FILES

This update package consists of 4 files in total (1 main update file + 3 auxiliary files).

1. Download the complete update package at this link:

   https://drive.google.com/drive/folders/1crZ-yMi-Fz1g08tgEawUbTdikatut1VZ?usp=drive_link

2. Copy the Main Update File (~178MB) to the root folder of your USB stick.
   - *Filename:* resmonpro-23.1.0a-for23.0.0-update.tar.bz2
3. Copy the remaining 3 Auxiliary Files to the root folder of the same USB stick.

RECOMMENDATIONS:

- It is advisable to download the files individually, as some browsers may rename them if downloaded all at once, which could cause the update to fail.
- If the file transfer to the USB fails for any reason (such as USB unplugging or I/O errors), it is recommended to restart the process—delete the already transferred files and reload them onto the USB drive.
- All files must be placed in the root folder of the USB drive.

INSTALL THE UPDATE

The update installation via the RESMON PRO FULL follows the same process as previous updates:

1. 1. Connect the USB stick containing the update files to the RESMON PRO FULL
2. 2. Log in as ADMIN, then navigate to the UPDATE section
3. 3. Select resmonpro-23.1.0a-for23.0.0-update.tar.bz2, and
4. 4. Click INSTALL.